



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
BINOCULAR INDIRECT OPHTHALMOSCOPE

Product Code :
SAMCOMEDICAL4004A044



Description :

BINOCULAR INDIRECT OPHTHALMOSCOPE Brilliant, Natural & Even illumination from a 6v 20w Halogen Lamp. Built in blue and red free filters & three step field of illumination. Compact and light weight for longer periods without fatigue. Precision optics and elegant mechanism, synchronizes the convergence and parallel adjustments. Also available in Fiber Optic Light system Contents: One condensing lens & carrying case, Teaching Mirror, S. Depressor, Fundus Chart, One spare Bulb. Head Band Indirect Ophthalmoscope with A494 Light Source & F.O. Cable C-450 Head Band Indirect Ophthalmoscope Halogen Complete Replacement Bulb OSRAM BINOCULAR INDIRECT OPHTHALMOSCOPE, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,