



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
OPHTHALMIC / ENT MICROSCOPE

Product Code :
SAMCOMEDICAL4004A048



Description :

OPHTHALMIC / ENT MICROSCOPE Magnification 5x to 25x continues. Eye piece pair SWF 15x. Inbuilt light source 15v/150w,twin bulb. Illumination through fibre optic cable 6mm 4ft. Maximum illumination 90,000lux. Fine focusing manual 54mm. Adjustable handle for handling. Working distance 200mm Heavy floor stand with five caster wheel. Optional Objective lens f-300mm, f-400mm. Beem splitter.C Mount with camera. Beem Splitter with Side observation tube. OPHTHALMIC / ENT MICROSCOPE, ENT Equipment exporters, ENT Equipment manufacturers



Samco Medical

Website: www.samcomedical.com,