



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ishihara Colour Vision Test

Product Code :
SAMCOMEDICAL240001A004



Description :

Tests for protan and deutan defects. Patients must identify numbers or patterns formed in a matrix of small coloured dots varying in colour and intensity (psudeo-isochromatic plates). Two editions available: 38 plate for detailed investigation - 24 plate for quick screening Ishihara Colour Vision Test, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,