



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Electric 7 Function

Product Code :
SAMCOMEDICAL1001A005

Description :

Samcomedical ICU Bed Electric 7 Function

Technical Specification :

Samcomedical ICU Bed Electric 7 Function with electrically operated back rest tilting 0-75, knee rest tilting 0-45, trendelenburg tilting 0-27, Reverse trendelenburg tilting 0-27, mattress base tilting to the left up to 35, tilting to the right up to 35. Size: L 2170mm W900mm

H 520-910mm(Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Separate Movable Nurse Control Panel can lock the buttons of Guard Rail Control panels and Foot Step control panel. High quality beautiful ABS engineering plastic head panel and foot panel equipped with bumper protection. Epoxy coated mild steel frame work and 4 section perforated top. Equipped with battery backup. Noiseless Castors with central locking system. Provision for I.V. Rod on both sides of the bed



Samco Medical

Website: www.samcomedical.com,