



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical ICU Bed Electric 7 Function With X-Ray  
Permeable Backrest

**Product Code :**  
SAMCOMEDICAL1001A006

**Description :**

Samcomedical ICU Bed Electric 7 Function With X-Ray Permeable Backrest

**Technical Specification :**

Samcomedical ICU Bed Electric 7 Function With X-Ray Permeable Backrest with electrically operated back rest tilting 0-80, knee rest tilting 0-35, trendelenburg tilting 0-20, Reverse trendelenburg tilting 0-20, mattress base tilting to the left up to 40, tilting to the right up to 40.

Size: L 2140mm

W 940mm

H 590-990mm (Adjustable Height). With X-ray permeable backrest. All functions controlled with DEWERT/T-motion motor from Germany/Taiwan Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Embedded controls on both sides of guard rails. Foot Step Control Panel under the bed frame. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. Equipped with battery backup. Noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

