



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

I.C.U. Multi-function Electric Bed With Weighing Scale

**Product Code :**

SAMCOMEDICAL1001A007

**Description :**

I.C.U. Multi-function Electric Bed With Weighing Scale

**Technical Specification :**

I.C.U. Multi-function Electric Bed With Weighing Scale Multi-Function ICU Bed electrically operated back rest tilting 0-65°, Foot rest tilting 0-43°, trendelenburg tilting 10°, Reverse trendelenburg tilting 10°, Height Adjustment (710mm to 530mm). Size: L 2330mm x W1050mm x H 530-710mm (Adjustable Height). All functions controlled with LINAK (Denmark) Power device (Actuator) In Build control panel on both side of Guard rail for Nurse/Patients to adjust bed position easily. Patient Lock-Out system for Nurse Staff at foot board. CPR position button for bed to reach 0° position quickly. Bed position data can be stored and recalled. Function of weighing scale, weight data can be stored and managed. Off-bed alarming. Removable High Quality ABS engineering Plastic head & foot panel. Size 1955mmX 900mm Four easy lifting Guard rail controlled by gas spring. Epoxy coated Mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 150mm dia swivel noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),