



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Bed Plain

**Product Code :**  
SAMCOMEDICAL1001A004

**Description :**

Hospital Bed Plain

**Technical Specification :**

Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work. Mounted on protective stumps. Bumper protection at 4 corners. Provision for I.V. Rod. Freight saving knock-down construction. Size: 2120 ½ 940 ½ 500 mm.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),