



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

I.C.U Bed Multifunction Electric Five Function

Product Code :

SAMCOMEDICAL1001A008

Description :

I.C.U Bed Multifunction Electric Five Function

Technical Specification :

I.C.U Bed Multifunction Electric Five Function Multi-Function ICU Bed with electrically operated back rest tilting 0-70°, Foot rest tilting 0-40°, trendelenburg tilting 12°, Reverse trendelenburg tilting 12°, Height Adjustment (710mm to 510mm). Size: L 2200mm x W1050mm x H 510-710mm (Adjustable Height). All functions controlled with LINAK (Denmark) power device (Actuator) In build control panel on both side of Guard rail for Nurse/Patients to adjust bed position easily. Patient Lock-Out system for Nurse Staff at foot board. CPR position button for bed to reach 0° position quickly. Removable High Quality ABS engineering Plastic head & foot panel. Sz 1945mm X 900mm Four easy lifting Guard rail controlled by gas spring. Epoxy coated Mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 150mm dia swivel noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed.



Samco Medical

Website: www.samcomedical.com,