



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Electric Five Function

Product Code :
SAMCOMEDICAL1001A009

Description :

Samcomedical ICU Bed Electric Five Function

Technical Specification :

Samcomedical ICU Bed Electric Five Function with electrically operated back rest tilting 0-80, knee rest tilting 0-35, trendelenburg tilting 0- 20, Reverse trendelenburg tilting 0-20

Size: L 2140mm

W 940mm

H 590-990mm (Adjustable Height). All functions controlled with DEWERT/T-motion motor from Germany/Taiwan Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. Equipped with battery backup. 125mm dia noiseless castors with Central braking system. Back rest & Knee rest can be adjusted simultaneously with single button control. Provision for I.V. Rod on both sides of the bed



Samco Medical

Website: www.samcomedical.com,