



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Electric 5 Function

Product Code :
SAMCOMEDICAL1001A010

Description :

Samcomedical ICU Bed Electric 5 Function

Technical Specification :

Samcomedical ICU Bed Electric 5 Function with electrically operated remote controlled back rest tilting 0-80, knee rest tilting 0-35, trendelenburg tilting 0-20, Reverse trendelenburg tilting 0-20. Size: L 2140mm

W 940mm

H 480-730mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. Collapsible side railings (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy). Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



Samco Medical

Website: www.samcomedical.com,