



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Mechanical Five Function

Product Code :
SAMCOMEDICAL1001A011

Description :

Samcomedical ICU Bed Mechanical Five Function

Technical Specification :

Samcomedical ICU Bed Mechanical Five Function with mechanically operated back rest, knee rest tilting 0-35, trendelenburg 0-20/ reverse trendelenburg 0-20 (TR/RTR) and height position. Size: L 2140mm

W 980mm

H 420-720mm (Adjustable Height). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. C Type design. Epoxy coated mild steel frame work and 4 section perforated top. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed.



Samco Medical

Website: www.samcomedical.com,