



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
I.C.U Bed, Hydraulic (Five Functions)

**Product Code :**  
SAMCOMEDICAL1001A012

**Description :**

I.C.U Bed, Hydraulic (Five Functions)

**Technical Specification :**

I.C.U Bed, Hydraulic (Five Functions) ICU bed Five-function with Hydraulic height adjustment by pressing UP/ DOWN pedal and mechanically operated back rest tilting 0-80° and Knee rest tilting 0-45°. Trendelenburg tilting 0-25°, Reverse trendelenburg tilting 0-25° is also control by switch at end foot. Size: L 2160mm x W 990mm x H 845-440mm (Adjustable Height). Removable High Quality ABS engineering Plastic head & foot panel Head panel and foot panel equipped with roller bumpers. 125mm dia noiseless castors with Central braking system



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),