



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical ICU Bed Hydraulic Four Function

**Product Code :**  
SAMCOMEDICAL1001A013

**Description :**

Samcomedical ICU Bed Hydraulic Four Function

**Technical Specification :**

Samcomedical ICU Bed Hydraulic Four Function with hydraulically operated back rest tilting 0-60,  
Height adjustable 520-810mm. Trendelenburg & Reverse Trendelenburg Position controlled by gear mechanic system tilting 0-11.  
Size: L 2240mm  
W 980mm  
H 520-810mm (Adjustable Height). Removable & interchangeable flat tube head panel and foot panel with . Epoxy coated mild steel frame work and 2 section Mesh Mattress Platform. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),