



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Electric Three Function

Product Code :
SAMCOMEDICAL1001A014

Description :

Samcomedical ICU Bed Electric Three Function

Technical Specification :

Samcomedical ICU Bed Electric Three Function with electrically operated remote controlled back rest 0-70, knee rest 0-35 and height position. Size: L 2160mm

W 940mm

H 520-770mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



Samco Medical

Website: www.samcomedical.com,