



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

I.C.U. Bed, Electric, Three Function

Product Code :

SAMCOMEDICAL1001A015

Description :

I.C.U. Bed, Electric, Three Function

Technical Specification :

I.C.U. Bed, Electric, Three Function Three-Function Electric Bed with electrically operated remote controlled back rest tilting 0-75°, Knee rest tilting 0-35°, Height Adjustment (640mm to 430mm). Size: L 2155mm x W 1050mm x H 430-640mm (Adjustable Height). All functions controlled with LINAK (Denmark) power device (Actuator) Removable High Quality ABS engineering Plastic head & foot panel, powder coated steel bed board with ventilated holes 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed



Samco Medical

Website: www.samcomedical.com,