



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical ICU Provision Bed Electric 3 Function

**Product Code :**  
SAMCOMEDICAL1001A016

**Description :**

Samcomedical ICU Provision Bed Electric 3 Function

**Technical Specification :**

Samcomedical ICU Provision Bed Electric 3 Function with electrically operated remote controlled back rest 0-70, knee rest 0-35 and height position. All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Size: L 2140mm

W 940mm

H 500-750mm (Adjustable Height). Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy). Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),