



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Electric Two Function

Product Code :
SAMCOMEDICAL1001A017

Description :

Samcomedical ICU Bed Electric Two Function

Technical Specification :

Samcomedical ICU Bed Electric Two Function with electrically operated remote controlled back rest 0-70, knee rest 0-35 position. All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Size: L 2140mm

W 940mm H 500mm. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors (all with individual brakes.) Provision for I.V. Rod on both sides of the bed



Samco Medical

Website: www.samcomedical.com,