



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Manual Three Function

Product Code :
SAMCOMEDICAL1001A018

Description :

Samcomedical ICU Bed Manual Three Function

Technical Specification :

Samcomedical ICU Bed Manual Three Function with manually operated back rest 0-80, knee rest 0-35 and height position.

Size: L 2160mm

W 940mm

H 530-810mm (Adjustable Height). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed



Samco Medical

Website: www.samcomedical.com,