



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical ICU Bed Provision Electric 5 Function

Product Code :
SAMCOMEDICAL1001A022

Description :

Samcomedical ICU Bed Provision Electric 5 Function

Technical Specification :

Samcomedical ICU Bed Provision Electric 5 Function with electrically operated remote controlled back rest tilting 0-80, knee rest tilting 0-35, trendelenburg tilting 0-20, Reverse trendelenburg tilting 0-20

Size: L 2160mm

W 940mm H 480-730mm (Adjustable Height). All functions controlled with Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with CPR device which allows nursing staff & doctors to simultaneously lower the back rest to a full flat position from any degree of elevation in order to save vital time in critical situation. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



Samco Medical

Website: www.samcomedical.com,