



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical ICU Bed Electric Weighing Scale

**Product Code :**  
SAMCOMEDICAL1001A023

**Description :**

Samcomedical ICU Bed Electric Five Function With Weighing Scale System

**Technical Specification :**

Samcomedical ICU Bed Electric Five Function With Weighing Scale System with electrically operated back rest tilting 0-80, knee rest tilting 0-35, trendelenburg tilting 0-20, Reverse trendelenburg tilting 0-20

Size: L 2160mm

W 940mm

H 540-770mm (Adjustable Height). All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with battery backup & with Weighing scale system. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Easy to operate built-in Control Panel on both sides of guard rails. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with bumpers protection. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),