



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical ICU Bed Electric Six Function

Product Code :

SAMCOMEDICAL1001A026

Description :

Samcomedical ICU Bed Electric Six Function With Sitting Position

Technical Specification :

Samcomedical ICU Bed Electric Six Function With Sitting Position with electrically operated sitting position 0-70, back rest tilting 0 - 75%, knee rest tilting 0-35, trendelenburg tilting 0 - 152, Reverse trendelenburg tilting 0 - 152

Size: 2155 x 950 x 570 / 780 mm (Adjustable Height). All functions controlled with Electric Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Head panel and foot panel equipped with safety lock and roller bumpers. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed



Samco Medical

Website: www.samcomedical.com,