



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Cold Boxs

Product Code :
SAMCOMEDICAL220002A001

Description :

Vaccine Storage Capacity : 23.3 Litres Weight Fully Loaded : 47.5 Kgs. Weight Empty (with empty ice pack) : 22.2 Kgs. Insulation Thickness : 100 mm Approx. External Dimensions(WxDxH) : 77x62x52 cms Approx. Internal Dimensions(WxDxH) : 52x37x28 cms Vaccine Storage Dimension : 45x30x18 cms Number of Ice Packs supplied : 50 Volume of Ice Packs : 0.30 Litres Conforming to WHO Specifications : E4/CB.2 Cold Life Requirement as per WHO at 43 $\frac{1}{2}$ C $\frac{1}{2}$ 1 $\frac{1}{2}$ C without opening : Minimum 120 hrs. Actual Cold Life without opening at 43 $\frac{1}{2}$ C $\frac{1}{2}$ 1 $\frac{1}{2}$ C : 141 hrs. 12 min. Cold Box, Vaccine Carrier Exporters, Vaccine Carrier Manufacturers



Samco Medical

Website: www.samcomedical.com,