



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Anti Bedsore Air Jet Type Mattress

Product Code :
SAMCOMEDICAL10011A009

Description :

Anti Bedsore Air Jet Type Mattress

Technical Specification :

Anti Bedsore Mattress Anti Bedsore Mattress (Air Jet Type). Tiny Holes on the surface for decreasing humidity. Replaces the foul air & reduce growth of microbes. Power supply: 220V/50Hz, 1A/250V, Air output: 3.5-4 liter/min. Pressure range: 50 mmHg 90 mmHg. Material: PVC Size:200cm(L) x 90cm(W) x 6.5cm (H) Low Noisy, Long life time benefit



Samco Medical

Website: www.samcomedical.com,