



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Bedside Screen Regular

Product Code :
SAMCOMEDICAL10016A006

Description :

Samcomedical Bedside Screen Regular

Technical Specification :

Bedside Screen - Regular Bedside Screens, 3 folds Size : 168H x 244W cms. Frame work made of ERW tubular pipe. Supplied with spring wire and curtains. Pre-treated and epoxy powder coated. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,