



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Saline Stand with IV Stand Economy

Product Code :

SAMCOMEDICAL10013A003

Description :

Samcomedical Saline Stand with IV Stand Economy

Technical Specification :

Saline Stand/IV Stand - Economy S.S. tubular pipe. Five legged heavy polypropylene base mounted on 5 cms castors. Double hook S.S. adjustable IV rod with a bakelite knob. Pre-treated and epoxy powder coated. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,