



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stand S.S.Tubular Pipe Economy

Product Code :
SAMCOMEDICAL10013A004

Description :

Saline Stand with IV Stand S.S.Tubular Pipe Economy

Technical Specification :

Saline Stand / IV Stand, Economy Saline Stand / IV Stand, Economy, S.S.tubular pipe. Five legged heavy polypropylene base mounted on 5 cm castors 4 hook adjustable S.S. I.V. Rod with a bakelite knob. Freight saving knock down construction.



Samco Medical

Website: www.samcomedical.com,