



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Delivery Table / Delivery Couch Stainless Steel

Product Code :

SAMCOMEDICAL10014A002

Description :

Delivery Table / Delivery Couch Stainless Steel

Technical Specification :

Delivery Table / Delivery Couch Stainless Steel Delivery Table / Bed, made in 3 sections. Size: 1680 $\frac{1}{2}$ 600 $\frac{1}{2}$ 800 mm. With adjustable Back rest & leg rest to varying positions. Framework made of Stainless Steel. Supplied with S.S. bowl and two leg holders. Mounted on protective stumps. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,