



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Obstetric Bed - Multi Function

Product Code :
SAMCOMEDICAL10014A004

Description :

Obstetric Bed - Multi Function

Technical Specification :

Obstetric Bed - Multi Function This table is designed for women's obstetric and gynaecological operation, Easy operation, safe and reliable, economical and quick. Operating positions - up / down, forward and backward tilt are controlled by smooth hand wheel. Comfortable for both doctor and Patient. Length and width for tabletop: 1800 x 600 mm. Height adjustment range: 650 - 895mm. Dimension of sections: backrest: 600 X 800 mm, Seat section: 600 x 400 mm, Leg section: 540 x 500 mm. Max. degree for backrest adjustment: more than 60 degree upwards, more than 10 degree downwards. Max. degree for tabletop forward and backward tilting: more than 8 degree(forward), more than 20 degree(backward).



Samco Medical

Website: www.samcomedical.com,