



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Obstetric Table - Multi Function

**Product Code :**  
SAMCOMEDICAL10014A006

**Description :**

Obstetric Table - Multi Function

**Technical Specification :**

Obstetric Table - Multi Function Obstetric table 3 section is designed for women's obstetric and gynecological examination, operation and delivery as well. It is Easy to operate, safe, reliable, and quick. Body section height is adjustable by hydraulic pump where as leg section height is manually adjusted. Back rest is adjusted by gas spring Size of the platform : Main table 1350ï¿½650mm, leg table 800ï¿½550mm. Height of the main table: ?650/850mm. height of the leg table: ?650/850mm.Backboard fold range: upward ?65ï¿½, downward ?5ï¿½. Platform tilt forward ?8ï¿½. Platform tilt backward ?20ï¿½ Both the main table and the leg table are mounted on castors for easy movement and can be locked up



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),