



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Obstetric Table Multi Function

Product Code :
SAMCOMEDICAL10014A007

Description :

Samcomedical Obstetric Table Multi Function

Technical Specification :

Obstetric Table - Multi Function Multi Function Obstetric table is designed for women's gynecological examination diagnosis, operation and delivery as well. Height adjustment, Trendelenburg and reverse trendelenburg of the whole table and the backrest are controlled by advanced linear motors through foot switch. Minimum height of the platform: 650mm. Height adjustment range: 240mm. Platform tilt forward 10°. Platform tilt backward 25°. Backboard fold range: upward 65°, downward 5°. Length of the main table: 1240mm. Width of the main table 600mm. Length of the leg table : 500mm. Width of the leg table : 600mm



Samco Medical

Website: www.samcomedical.com,