



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Labour & Delivery Bed (With Accessories)

**Product Code :**  
SAMCOMEDICAL10014A008

**Description :**

Labour & Delivery Bed (With Accessories)

**Technical Specification :**

Labour & Delivery Bed (With Accessories) Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes. Leg section can be completely recessed under the body section. Size : Body section: 1000(L) x 800(W) x 750(H)mm., Leg section: 900(L) x 800(W) x 750(H)mm. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Mattress-pad 2 sections. Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width. Epoxy coated tubular steel frame. Mattress-pad density : 21-25 kg/m<sup>3</sup>. Carrying capacity: 150 - 180 kg. Mattress-pad: Approx. 60(H)mm. Supplied with S.S. bowl. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),