



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Gynecology Table - DISMANTLED(With Accessories)

Product Code :
SAMCOMEDICAL10014A009

Description :

Gynecology Table DISMANTLED(With Accessories)

Technical Specification :

Gynecology Table - DISMANTLED(With Accessories) Made in 3 sections, adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips. Size : Table 1800L x 600W x 800H mm. Tubes : 30 x 1.5 mm. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Support ring at proximity of seat section to support round basin-bowl. Supplied with 2 leg holders with canvas straps adjustable in height and width. Epoxy coated tubular steel frame. Tubes : 30 x 1.5 mm. Mattress-pad : density : 21 - 25 kg/m³ ± ½. Carrying capacity : 150-180 kg. Mattress-pad : 30H mm. Knockdown construction. Supplied with 4 liter S.S. Bowl. Note : Available with knee crutches on request.



Samco Medical

Website: www.samcomedical.com,