



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Delivery Bed

Product Code :
SAMCOMEDICAL10014A010

Description :

Delivery Bed

Technical Specification :

Delivery Bed Size : 1800L x 550W x 800H mm. Adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Supplied with 2 knee crutches, adjustable in height and width. Framework made of mild steel epoxy coated finish with upholstered top. Supplied with S.S. bowl. Knockdown construction



Samco Medical

Website: www.samcomedical.com,