



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Labour and Delivery Bed

**Product Code :**  
SAMCOMEDICAL10014A011



**Description :**

Samcomedical Labour and Delivery Bed With Accessories

**Technical Specification :**

Labour & Delivery Bed (With Accessories) Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes. Size : Body section: 1000L x 800W x 750H mm., Leg section: 900L x 800W x 750H mm. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width. Epoxy coated tubular steel frame. Carrying capacity: 150 - 180 kg. Supplied with S.S. bowl. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),