



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Delivery Table Telescopic (Fixed Height)

Product Code :
SAMCOMEDICAL10014A012

Description :

Delivery Table Telescopic (Fixed Height)

Technical Specification :

Delivery Table Telescopic (Fixed Height) Size : 182L x 68W x 76H cms. Frame work made of MS Tube. Pretreated & Epoxy powder coated finish 3 Section Top made of Stainless Steel Sheet. Legs mounted on PVC stumps Trendlenburg positions adjustable by pneumatic gas spring system. Backrest section adjustable on ratchet. Leg end section can slide under the main section Height adjustable pair of foam padded knee crutches. Removable Safety Side railing provided on 3 side of Table. Supplied with waste collecting Stainless Steel Bowl on S.S. Bracket and one Stainless Steel I V Rod.



Samco Medical

Website: www.samcomedical.com,