



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Delivery Table Telescopic

Product Code :

SAMCOMEDICAL10014A014

Description :

Samcomedical Delivery Table Telescopic Adjustable Height

Technical Specification :

Delivery Table Telescopic (Adjustable Height) Size : 182L x 68W x 60-80H cms. Tubular frame work mounted on swiveling castors. Pre treated & Epoxy powder coated finish. 3 Section Top made of Stainless Steel Sheet. Trendlenburg positions adjustable by pneumatic gas spring system. Hydraulically Height adjustable by foot paddle from side. Backrest section adjustable on ratchet. Leg end section can slide under the main section Height adjustable pair of foam padded knee crutches. Removable Safety Side railing provided on 3 side of Table. Supplied with waste collecting Stainless Steel Bowl on S.S. Bracket and one Stainless Steel I V Rod.



Samco Medical

Website: www.samcomedical.com,