



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Labour Table / Delivery Bed

Product Code :
SAMCOMEDICAL10014A015

Description :

Labour Table / Delivery Bed

Technical Specification :

Labour Table / Delivery Bed Made in 3 sections with independent controls and trendelburg position. Adjustable back rest and leg section with mechanism to varying positions, capable of tilting for obstetric and gynecological use. Size of assembled table : Approx. 1800L x 550W x 800H mm. Framework made of chrome plated rectangular steel tube with upholstered top. Supplied completely with S.S. bowl and two standard padded leg holders. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,