



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Medicine Trolley

Product Code :
SAMCOMEDICAL10015A010

Description :

Medicine Trolley

Technical Specification :

Medicine Trolley Two side Trolley designed to supply medicines. (Front & back same) Size: 850 ½ 520 ½ 1045 mm. 2 drawers with dividers suitably designed for keeping medicines. 3 Multi purpose drawers. Equipped with waste bin, needle disposable container, file cassette & guard rails. Equipped with lock key system. 4 swivel noiseless castors with brakes.



Samco Medical

Website: www.samcomedical.com,