



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Medicine Trolley

Product Code :
SAMCOMEDICAL10015A011

Description :

Samcomedical Medicine Trolley

Technical Specification :

Medicine Trolley Trolley made of Alloy, ABS plastic and stainless steel Size : 850mm x 520mm x 1010mm 6 drawers, can be pulled out from both sides, self locking when closed. 50 medicine boxes in all drawers History card on each medicine box. With waste bins-2 pcs, Transparent file box-1 pcs, sliding side board, and five hooks hanging device mechanism. 4 swivel noiseless castor of 5" (2 with brakes).



Samco Medical

Website: www.samcomedical.com,