



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Emergency Trolley

Product Code :
SAMCOMEDICAL10015A012

Description :

Samcomedical Emergency Trolley

Technical Specification :

Emergency Trolley Trolley designed to manage patient in emergency. Size: 850 520 1045 mm. 5 drawers made of ABS with dividers. Equipped with waste bin, needle disposable container, multi bin container, file cassette, guard rails, Sliding shelf & defibrillator shelf. Equipped with Oxygen Cylinder Holder, I.V. Pole & CPR Board. Equipped with lock key system. (Central locking mechanism) 4 swivel noiseless castors of 5" (2 with brakes).



Samco Medical

Website: www.samcomedical.com,