



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Treatment Trolley

**Product Code :**  
SAMCOMEDICAL10015A015

**Description :**

Treatment Trolley

**Technical Specification :**

Treatment Trolley 90% made of ABS. Size: 818 x 525 x 985 mm. 2 concealed working tops under table top. 2 drawers & 2 cupboards. With medicine containers/drawers. With 2 waste bins & I.V. Pole. 4 swivel castors (2 with brakes). ABS table top , Sz. 700 x490mm



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),