



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Emergency Two Side Trolley

Product Code :
SAMCOMEDICAL10015A018

Description :

Emergency Two Side Trolley

Technical Specification :

Emergency Trolley Two side Trolley designed to supply medicines. (Front & back same) Size: 850 $\frac{1}{2}$ 515 $\frac{1}{2}$ 985 mm. 4 drawers with dividers suitably designed for keeping medicines. 1 Multi purpose drawers. Equipped with IV pole, waste bin, needle disposable container, file cassette & guard rails. Equipped with lock key system. 4 swivel castors (2 with brakes).



Samco Medical

Website: www.samcomedical.com,