



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Medicine Trolley, Superior S.S.

Product Code :
SAMCOMEDICAL10015A024

Description :

Medicine Trolley, Superior S.S.

Technical Specification :

Medicine Trolley, Superior S.S. Trolley made of stainless steel. Size: 640 ½ 450 ½ 850 mm. Top surface with 3 side guard railings. 2 S.S. lockable drawers. 1 S.S. Shelf below the drawers. 4 swivel castors (2 with brakes)



Samco Medical

Website: www.samcomedical.com,