



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dressing Trolley, ABS

Product Code :
SAMCOMEDICAL10015A025

Description :

Dressing Trolley, ABS

Technical Specification :

Dressing Trolley, ABS Size: 820 ½ 510 ½ 970 mm. 2 ABS shelves with guard rails on 3 sides. Main frame work made of Stainless Steel. With 1 waste bin & 1 S.S. Bowl. 1 drawer below the first shelf. 4 swivel castors (2 with brakes).



Samco Medical

Website: www.samcomedical.com,