



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Dressing Trolley ABS

Product Code :
SAMCOMEDICAL10015A026

Description :

Samcomedical Dressing Trolley ABS

Technical Specification :

Dressing Trolley, ABS Size: 820 510 970 mm. 2 ABS shelves with guard rails on 3 sides. With dressing tray attached at the top and epoxy coated frame work. With 1 waste bin & 1 S.S. Bowl. 1 drawer below the first shelf. 4 swivel castors (2 with brakes).



Samco Medical

Website: www.samcomedical.com,