



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dressing Trolley, S.S.

Product Code :
SAMCOMEDICAL10015A027

Description :

Dressing Trolley, S.S.

Technical Specification :

Dressing Trolley, S.S. Trolley made of stainless steel. Size: 660 ½ 440 ½ 925 mm. 2 S.S. shelves with guard rails on 3 sides. With 1 waste bin & 1 S.S. Bowl. S.S. Dressing tray attached at the top. 1 S.S. drawer below the first shelf. 4 swivel castors (2 with brakes).



Samco Medical

Website: www.samcomedical.com,