



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Trolley For Records ABS

Product Code :

SAMCOMEDICAL10015A034

Description :

Samcomedical Trolley For Records ABS

Technical Specification :

Trolley For Records, ABS Trolley framework made up of Epoxy coated steel with 4 Aluminium column. ABS Working surface & base. Suitable for 25 clips. Size 350x400x1025mm. Equipped with Lock & Key. 4 Swivel castors with brakes.



Samco Medical

Website: www.samcomedical.com,