



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Mayo Type Trolley

Product Code :
SAMCOMEDICAL10015A040

Description :

Mayo Type Trolley

Technical Specification :

Mayo Type Trolley Rectangular CRC tubes mounted on 5cms castors, height adjustable by side knob. Height adjustable approx. 800 to 1200H mm. Supplied with Mayo tray, stainless steel. Size: 560L x 400W x 20H mm. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,