



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Dressing Trolley

Product Code :
SAMCOMEDICAL10015A043

Description :

Samcomedical Dressing Trolley

Technical Specification :

Dressing Trolley Size : 76L x 51W x 82H cms. Frame work made of ERW tubular pipe of 31mm dia. mounted on four 10 cms. dia castors. Two fixed type S.S. Heavy gauge shelves. Supplied complete with bucket and basin. Pre-treated and epoxy powder coated. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,