



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Dressing Trolley - Stainless steel

**Product Code :**  
SAMCOMEDICAL10015A044

**Description :**

Dressing Trolley - Stainless steel

**Technical Specification :**

Dressing Trolley - Stainless steel Size : 76L x 51W x 82H cms. Frame work made of Stainless steel mounted on four 10 cms. dia castors. Two fixed type S.S. heavy gauge shelves. Completely supplied with bucket and basin. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),